

# NAUTIAGENE

Product Information

## NautiaZ Bacteria/Fungi RNA Mini Kit

(100/300 prep)

Cat.No. : NGBRZ-S100/NGBRZ-S300

Sample : 10<sup>9</sup> Bacteria cells

10<sup>8</sup> Fungi cells

Yield : Up to 30 µg

NGCRZ-S100	NautiaZ Cell/Blood Total RNA Mini Kit (100 prep)
NGCRZ-S300	NautiaZ Cell/Blood Total RNA Mini Kit (300 prep)
NGBRZ-S100	NautiaZ Bacteria/Fungi RNA Mini Kit (100 prep)
NGBRZ-S300	NautiaZ Bacteria/Fungi RNA Mini Kit (300 prep)
NGTRZ-S100	NautiaZ Tissue Total RNA Mini Kit (100 prep)
NGTRZ-S300	NautiaZ Tissue Total RNA Mini Kit (300 prep)
NGPRZ-S100	NautiaZ Plant Total RNA Mini Kit (100 prep)
NGPRZ-S300	NautiaZ Plant Total RNA Mini Kit (300 prep)
NGMRZ-S050	NautiaZ microRNA Mini Kit (50 prep)
NGVN-S100	Nautia Viral Nucleic Acid Extraction Kit (100 prep)
NGVN-S300	Nautia Viral Nucleic Acid Extraction Kit (300 prep)

Ver. 2020-01

## Contents

	NGBRZ-S100T	NGBRZ-S100	NGBRZ-S300
BRO Buffer	1 ml	25 ml	65 ml
BRA Buffer	2 ml	45 ml	125 ml
W1 Buffer	2 ml	45 ml	125 ml
W2 Buffer*	300 µl x2	15 ml	25 ml x2
Elution Buffer	1 ml	10 ml	30 ml
RZ Column	4 pcs	100 pcs	300 pcs
Collection Tube	4 pcs	100 pcs	300 pcs
User Manual	1	1	1

\*Add 1.2 ml x2 / 60 ml / 100 ml x2 ethanol (96-100%) to W2 Buffer prior to the initial use

## Buffer Preparation

- Add ethanol (96-100%) to the Wash Solution prior to first use:

	NGBRZ-S100T	NGBRZ-S100	NGBRZ-S300
W2 Buffer ethanol (96 ~ 100%)	300 µl x2 1.2 ml x2	15 ml 60 ml	25 ml x2 100 ml x2

## Additional Requirements


1. β-Mercaptoethanol
  2. RNase-free microcentrifuge tubes
  3. 70% ethanol
- For Bacteria sample:**
4. Lysozyme Buffer (20 mg/ml lysozyme; 20 mM Tris-HCl; 2 mM EDTA; 1% TritonX-100; pH 8.0, prepare the lysozyme buffer immediately prior to use)
- For Fungus sample:**
4. Lyticase or Zymolase
  5. Sorbitol buffer (1.2 M sorbitol; 10 mM CaCl<sub>2</sub>; 0.1M Tris-Cl pH 7.5; 35 mM mercaptoethanol)

## Important Notes

1. Buffer contains chaotropic salt is harmful and irritant agent.
2. Use sterile, RNase-free pipet tips and microcentrifuge tube. Wear a lab coat and disposable gloves to prevent RNase contamination.
3. Make sure the starting sample amount is under the limit.
4. Add ethanol (96-100%) to W2 Buffer prior to the initial use.
5. All centrifugations should be carried out in a table-top microcentrifuge at >12000 × g (~14000 × g).

## For Fungus

STEP	PROCEDURE
1 Sample prepare	<b>Transfer up to 10<sup>8</sup> fungus cells</b> to a microcentrifuge tube ( <i>not provided</i> ). Centrifuge at 6,000 x g for 5 mins. Remove the supernatant completely. Resuspend the cell with <b>600 µl Sorbitol Buffer</b> by pipetting. <b>Add 200 U of lyticase or zymolase.</b> Incubate at 30°C for 30 minutes.
	Centrifuge the mixture at 2,000 x g for 10 mins. Remove the supernatant completely.
	<b>Add 200 µl BRO Buffer</b> to the tube and resuspend the cell pellet by vortexing or pipetting. Incubate at room temperature for 5 minutes.
	<b>Add 300 µl BRA Buffer</b> and <b>3 µl β-Mercaptoethanol</b> , mix by vortexing. Incubate at room temperature for 5 minutes.



## For Gram- Bacteria

STEP	PROCEDURE
1 Sample prepare	Transfer up to $10^9$ bacteria cells to a microcentrifuge tube ( <i>not provided</i> ). Centrifuge at 12,000 x g for 1 min and remove the supernatant completely.
	<b>Add 200 <math>\mu</math>l BRO Buffer</b> to the tube and resuspend the cell pellet by vortexing or pipetting. Incubate at room temperature for 5 minutes.
	<b>Add 300 <math>\mu</math>l BRA Buffer</b> and <b>3 <math>\mu</math>l <math>\beta</math>-Mercaptoethanol</b> , mix by vortexing. Incubate at room temperature for 5 minutes.

## For Gram+ Bacteria

STEP	PROCEDURE
1 Sample prepare	Transfer up to $10^9$ bacteria cells to a microcentrifuge tube ( <i>not provided</i> ). Centrifuge at 12,000 x g for 1 min and remove the supernatant completely.
	<b>Add 200 <math>\mu</math>l Lysozyme Buffer</b> (20mg/ml lysozyme) to the tube and resuspend the cell pellet by vortexing or pipetting. Incubate at room temperature for 10 minutes.
	<b>Add 300 <math>\mu</math>l BRA Buffer</b> and <b>3 <math>\mu</math>l <math>\beta</math>-Mercaptoethanol</b> , mix by vortexing. Incubate at room temperature for 5 minutes.

## PURIFICATION PROTOCOLS

STEP	PROCEDURE
2 Cell Lysis	Centrifuge at 14,000 x g for 10 minutes. Transfer the supernatant to a clean microcentrifuge tube.
3 RNA Binding	<b>Add 500 <math>\mu</math>l 70% ethanol</b> and shake vigorously. Place a <b>RZ column</b> in Collection Tube. Transfer the sample mixture (up to 700 $\mu$ l once) to RZ column and centrifuge 1 minute at 14,000 x g. Discard the flow-through and place RZ Column back in the Collection tube.
4-1 Wash	<b>Add 400 <math>\mu</math>l W1 Buffer</b> to RZ Column. Centrifuge at 14,000 x g for 30 seconds. Discard the flow-through and place <b>RZ Column</b> back in the Collection tube.
Optional Step: DNase	<i>If DNA-free RNA is required, perform this optional step.</i> <b>Add 150 <math>\mu</math>l W2 Buffer (ethanol added)</b> into the RZ column. Centrifuge at full speed (16,000 x g) for 30 seconds. Discard the flow-through and place the RZ Column back in the Collection Tube.
	For each isolation reaction , premix <b>80 <math>\mu</math>l DNase I Incubation Buffer</b> with <b>2 <math>\mu</math>l DNase I</b> in a new sterile tube ( <b>Do not vortex!</b> ). <b>Add 82 <math>\mu</math>l of the DNase I solution</b> into the center of the RZ Column membrane and incubate at room temperature for 15 min.

4-2 Wash	<b>Add 600 <math>\mu</math>l W2 Buffer</b> (ethanol added) to <b>RZ Column</b> . Centrifuge at 14,000 x g for 30 seconds. Discard the flow-through and place <b>RZ Column</b> back in the Collection tube.
5 Dry	Centrifuge at 14,000 x g for 2 minutes to dry the column.
6 Elution	Place <b>RZ Column</b> to a clean 1.5 ml microcentrifuge tube ( <i>not provided</i> ). <b>Add 50-100 <math>\mu</math>l of preheated Elution Buffer (75°C)</b> into the center of the column matrix.  Stand at room temperature for 2 minutes. Centrifuge at 14,000 x g for 2 minutes to elute purified RNA.
7 Pure RNA	Store the RNA fragment at -80 °C.