

NAUTIAGENE

Product Information

NautiaZ microRNA Mini Kit

(50 prep)

Cat.No. : NGMRZ-S050
Sample : 10⁷ Culture cells
300 µl blood
10⁹ Bacteria cells
10⁸ Fungi cells
Yield : Up to 30 µg



NGCRZ-S100	NautiaZ Cell/Blood Total RNA Mini Kit (100 prep)
NGCRZ-S300	NautiaZ Cell/Blood Total RNA Mini Kit (300 prep)
NGBRZ-S100	NautiaZ Bacteria/Fungi RNA Mini Kit (100 prep)
NGBRZ-S300	NautiaZ Bacteria/Fungi RNA Mini Kit (300 prep)
NGTRZ-S100	NautiaZ Tissue Total RNA Mini Kit (100 prep)
NGTRZ-S300	NautiaZ Tissue Total RNA Mini Kit (300 prep)
NGPRZ-S100	NautiaZ Plant Total RNA Mini Kit (100 prep)
NGPRZ-S300	NautiaZ Plant Total RNA Mini Kit (300 prep)
NGMRZ-S050	NautiaZ microRNA Mini Kit (50 prep)
NGVN-S100	NautiaZ Viral Nucleic Acid Extraction Kit (100 prep)
NGVN-S300	NautiaZ Viral Nucleic Acid Extraction Kit (300 prep)

Contents

	NGMRZ-S050T	NGMRZ-S050
ML Buffer	2 ml x2	60 ml
MO Buffer	1 ml	15 ml
MA Buffer	2 ml	25 ml
W1 Buffer	2 ml	25 ml
W2 Buffer*	300 ul x2	15 ml
Elution Buffer	1 ml	10 ml
Red Column	4 pcs	50 pcs
Blue Column	4 pcs	50 pcs
Collection Tube	8 pcs	100 pcs
User Manual	1	1

*Add 1.2 ml x2 / 60 ml ethanol (96-100%) to W2 Buffer prior to the initial use

Buffer Preparation

- Add ethanol (96-100%) to the Wash Solution prior to first use:

	NGMRZ-S050T	NGMRZ-S050
W2 Buffer ethanol (96 ~ 100%)	300 ul x2 1.2 ml x2	15 ml 60 ml

Optional Step (DNA Residue Degradation)

Add 5 µl of 10X Reaction Buffer and 5 µl of DNase I, Amplification Grade, and incubating to the 50 µl final product. Let stand for 15 minutes at room temperature. Inactivate the DNase by adding 5 µl of the Stop Solution and heating at 70 °C for 10 minutes. (Suggesting product: Sigma #AMPD1-1KT, DNase I, Amplification Grade)

Description

NautiaZ microRNA Mini Kit provides a fast, simple, and cost-effective method for isolation of microRNA from whole blood, mammalian cells and bacterial cells. Detergents and chaotropic salt are used to lyse cells and inactivate RNase. The specialized high-salt buffering system further allows microRNA bases to bind to the glass fiber matrix of the spin column while contaminants pass through the column. Impurities are efficiently washed away, and pure microRNA is eluted with Elution Buffer without phenol extraction or alcohol precipitation needs. microRNA purified with NautiaZ microRNA Mini Kit is suitable for a variety of routine applications including gene regulation and functional The entire procedure can be completed within 25-40 minutes.

Important Notes

1. Buffer contains chaotropic salt is harmful and irritant agent.
2. Use sterile, RNase-free pipet tips and microcentrifuge tube.
Wear a lab coat and disposable gloves to prevent RNase contamination.
3. Make sure the starting sample amount is under the limit.
4. Add ethanol (96-100%) to W2 Buffer prior to the initial use.
5. All centrifugations should be carried out in a table-top microcentrifuge at >12000 × g (~14000 × g).

Additional Requirements

1. β-Mercaptoethanol
2. RNase-free microcentrifuge tubes
3. 70% ethanol

For Bacteria sample:

1. Lysozyme Buffer (20 mg/ml lysozyme; 20 mM Tris-HCl; 2 mM EDTA; 1% TritonX-100; pH 8.0, prepare the lysozyme buffer immediately prior to use)

For Fungus sample:

1. Lyticase or Zymolase
2. Sorbitol buffer (1.2 M sorbitol; 10 mM CaCl₂; 0.1M Tris-Cl pH 7.5; 35 mM mercaptoethanol)

Prepare



Cell Lysis



Binding



Wash



Dry



Elution

For Culture Cell

STEP	PROCEDURE
1 Sample prepare	Transfer up to 10^7 Culture mammalian cells to a microcentrifuge tube (<i>not provided</i>). Centrifuge at $6,000 \times g$ for 1 min and remove the supernatant completely.
	Add 100 μl ML Buffer to the tube and resuspend the cell pellet by vortexing or pipetting.

For Blood

STEP	PROCEDURE
1 Sample prepare	Collect blood in EDTA-Na2 treated collection tubes. Transfer up to 300 μ l blood to a 1.5 ml tube. Add 900 μl ML Buffer , then mix by insertion. Incubate the mixture on ice for 10 minutes, and invert every 5 minutes.
	Centrifuge at $4,000 \times g$ for 5 minutes at 4°C . Remove the supernatant completely. Add 100 μl ML Buffer and resuspend the pellet by pipetting.

For Gram+ Bacteria

STEP	PROCEDURE
1 Sample prepare	Transfer up to 10^9 bacteria cells to a microcentrifuge tube (<i>not provided</i>). Centrifuge at $12,000 \times g$ for 1 min and remove the supernatant completely.
	Add 200 μl Lysozyme Buffer (20mg/ml lysozyme) to the tube and resuspend the cell pellet by vortexing or pipetting. Incubate at room temperature for 10 minutes.

For Gram- Bacteria

STEP	PROCEDURE
1 Sample prepare	Transfer up to 10^9 bacteria cells to a microcentrifuge tube (<i>not provided</i>). Centrifuge at $12,000 \times g$ for 1 min and remove the supernatant completely.
	Add 100 μl MO Buffer to the tube and resuspend the cell pellet by vortexing or pipetting. Incubate at room temperature for 5 minutes.

For Fungus

STEP	PROCEDURE
1 Sample prepare	Transfer up to 10^8 fungus cells to a microcentrifuge tube. Centrifuge at $6,000 \times g$ for 5 mins. Remove the supernatant completely. Resuspend the cell with 600 μl Sorbitol Buffer by pipetting. Add 200 U of lyticase or zymolase . Incubate at 30°C for 30 minutes.
	Centrifuge the mixture at $2,000 \times g$ for 10 mins. Remove the supernatant completely.
	Add 100 μl MO Buffer to the tube and resuspend the cell pellet by vortexing or pipetting. Incubate at room temperature for 5 minutes.

PURIFICATION PROTOCOLS

STEP	PROCEDURE
2 Cell Lysis	Add 300 μl MA Buffer and 3 μl β-Mercaptoethanol , mix by vortexing. Incubate at room temperature for 5 minutes.

	Centrifuge at $14,000 \times g$ for 10 minutes. Transfer the supernatant to a clean microcentrifuge tube.
3 Protein remove	Add 300 μl 70% ethanol and mix by vortexing. Place a Red column in Collection Tube. Transfer the sample mixture (up to 700 μ l once) to Red column and centrifuge 1 minute at $14,000 \times g$. Transfer the the flow-through to a new 1.5 ml microcentrifuge tube. Add 600 μl 70% ethanol and mix by vortexing. Place a Blue column in Collection Tube.
4 RNA Binding	Transfer the sample mixture (up to 700 μ l once) to Blue column and centrifuge 1 minute at $14,000 \times g$. Discard the flow-through and place Blue Column back in the Collection tube.
5-1 Wash	Add 400 μl W1 Buffer to Blue Column. Centrifuge at $14,000 \times g$ for 30 seconds. Discard the flow-through and place Blue Column back in the Collection tube.
5-2 Wash	Add 600 μl W2 Buffer (ethanol added) to RZ Column . Centrifuge at $14,000 \times g$ for 30 seconds. Discard the flow-through and place RZ Column back in the Collection tube.
6 Dry	Centrifuge at $14,000 \times g$ for 2 minutes to dry the column.
7 Elution	Place Blue Column to a RNase-free 1.5 ml tube. Add 50 μl of preheated Elution Buffer (75°C) into the center of the column matrix.
	Stand at room temperature for 2 minutes. Centrifuge at $14,000 \times g$ for 2 minutes to elute RNA.
8 Pure RNA	Store the RNA fragment at -80°C .